



qLabs® ElectroMeter (Q3)

for Professional Point of Care



MICROPOINT

To find out more, please email: info@operonbiotech.com






qLabs® ElectroMeter (Q3)

for Professional Point of Care

REF Q-3 Plus



Diagnostics on Demand for POC

Touchscreen handheld rapid testing platform with option of eStation II docking system

Features

- Reliable laboratory results in 30 -100 seconds for PT and 4 - 7 minutes for aPTT
- 10 µL of fingertip whole blood, fresh venous whole blood or citrate whole blood
- On board quality control
- External liquid QC (LQC) available
- Built-in temperature controlled testing
- Built-in QR code scanner
- No calibration or instrument manipulation required

Current Test Menu

- PT-INR
- aPTT
- PT/aPTT Combo

ABOUT MICROPOINT

Micropoint Biotechnologies Co., Ltd. is a global provider of Point of Care Testing (POCT) systems and solutions. We are improving patient care by providing rapid, accurate, reliable, and low cost diagnostic products. Micropoint's patented microfluidic technologies and three proprietary POCT platforms enable testing at the same high quality level as a central lab, but with improved speed and convenience. Micropoint's products are manufactured in cGMP and ISO 13485 certified facilities.

MFD By:

MICROPOINT

Micropoint Biotechnologies, Inc.
3521 Leonard Court
Santa Clara, CA 95054
United States

Interface

Display

3.5" 320×480 Color LCD

Barcode scanner

QR code scanner

Software

Language

Multi-language menu

English, Czech, German, French, Italian,

Spanish, Portuguese, Dutch, Slovak, Romanian,

Norwegian, Swedish, Danish, Polish, Greek

Memory Capacity

Patient records

Stores up to 2,000 patient records
and 500 raw data

Performance

Measurable Range

INR range 0.5 to 8.0

aPTT from 20 sec to 420 sec

Specifications

Dimensions

148mm x 70mm x 26mm (L x W x H)

Weight

203 g

Operating conditions

Temperature: 10 °C - 35 °C

Relative humidity: 10 % - 90 %

Power supply

100V-240V AC

Built-in 2,500mAh lithium polymer battery

Imp & Mkt in India by

📍 **Operon Biotech & Healthcare**
#37, A-1, B-1, B-2, 'SnehaJyothi',
Srinivagilu Main Road
Koramangala Ring Road
Bengaluru - 560 047.

☎ 080-42027704 | Fax: 25701522

✉ operon@operonbiotech.com

www.operonbiotech.com

